# Geratherm<sup>®</sup>

## *active control*

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CE 0197

Digital wrist blood pressure monitor

# Geratherm<sup>®</sup>

*active control*

EN

**INSTRUCTION FOR USE**  
**GT-1215**

CE 0197

## Contents

|                                                              |    |
|--------------------------------------------------------------|----|
| <b>Indications for use</b>                                   | 37 |
| <b>Contraindications</b>                                     | 37 |
| <b>Precautions</b>                                           | 38 |
| <b>Cleaning</b>                                              | 39 |
| <b>Maintenance and Disposal</b>                              | 39 |
| <b>Storage</b>                                               | 40 |
| <b>Warranty</b>                                              | 40 |
| <b>Information You Should Know Before Operating the Unit</b> | 40 |
| <b>Measurement Principle</b>                                 | 43 |
| <b>About the Unit</b>                                        |    |
| Information about the Device                                 | 44 |
| Explanation of Display                                       | 45 |
| <b>Preparation for Use</b>                                   |    |
| Power Supply and Charging Power                              | 46 |
| Activate your Blood Pressure Monitor                         | 49 |
| Setting Date and Time                                        | 49 |
| Attaching the Cuff                                           | 51 |
| Measuring Posture                                            | 51 |
| <b>How to Operate the Unit</b>                               |    |
| Measuring Blood Pressure                                     | 52 |
| Retrieving Readings                                          | 54 |
| Deleting Readings                                            | 55 |
| <b>Care and Maintenance</b>                                  | 56 |
| <b>Error Messages</b>                                        | 58 |
| <b>Specifications</b>                                        | 59 |
| <b>Quality Standard</b>                                      | 60 |
| <b>Symbol Index</b>                                          | 60 |
| <b>Information on electromagnetic compatibility (EMC)</b>    | 61 |
| <b>Annex</b>                                                 | 62 |

Thank you for choosing the Geratherm® active control blood pressure monitor.

This unit uses an oscillometric measurement method in order to measure systolic and diastolic blood pressure, as well as the heart rate. The measurement is conducted on the wrist and all values appear on an LCD display screen.

**Before using your device for the first time, please read this instruction manual carefully.**

## Indications for use

This device is intended for non-invasive measuring and monitoring of arterial blood pressure. This device should only be used by adults over 18 years of age with a wrist circumference between 13.5 and 21.5 cm.

This unit was developed for home use exclusively.

EN

## Contraindications

This device is not intended for infants, pregnant women, patients with implanted electronic devices, patients with pre-eclampsia, premature ventricular beats, atrial fibrillation, peripheral artery disease and patients who are undergoing intravascular therapy or who have an arteriovenous shunt, or persons who have had a mastectomy. A doctor must be consulted prior to using the device on any of the persons listed above, and the suitability of the measuring method must be assessed.

This device is not suitable for continuous monitoring during medical emergencies and surgery. This device cannot be used simultaneously with HR surgical devices. The device is not intended for the transportation of patients outside medical establishments.

It is not intended for use on extremities other than the wrist or for functions other than obtaining a blood pressure measurement.

Do not wrap the cuff around an arm in which an arteriovenous shunt has been implanted or to which other medical devices are attached, as this could lead to temporary functional failure of the medical devices being used.

If you are allergic to polyester, nylon or plastic, please do not use this device.

The manufacturer cannot be held liable for any damage caused by improper use or failure to comply with the safety instructions in the instruction manual.

## Precautions

EN

- The device is intended for private use only.
- These instructions and the device are no substitute for visiting your doctor. Neither the information contained here nor this device may be used to diagnose or treat health problems or to prescribe medicines. If you have a medical problem or suspect that you may have such a problem, please seek immediate advice from your doctor.
- Do not mistake self-monitoring for self-diagnosis. This device enables you to monitor your blood pressure. Do not undertake any therapeutic measures on the basis of measurements you have taken yourself.
- Before using the device, ensure that it functions safely and is in perfect condition. Check the device. Do not use the device if it is damaged in any way. Continuing to use a damaged device can lead to injuries, false readings and serious danger.
- Keep the unit out of the reach of toddlers, children and pets. The inhalation or swallowing of small parts is dangerous and could be fatal.
- Please use only accessories that have been approved by the manufacturer. Otherwise, damage may be caused to the unit, the user may suffer injury or inaccurate measurements may occur.
- Please operate the unit under the environmental conditions described in the instructions. Otherwise, the performance and service life the blood pressure monitor will be adversely affected.
- This blood pressure monitor is not an AP/APG device and is not suitable for use in the presence of inflammable anaesthetic gases with air containing oxygen or nitrous oxide.
- If the cuff pressure reaches 300 mmHg (40 kPa), the air in the cuff is automatically released. If the cuff does not empty when the pressure reaches 300 mmHg (40 kPa), remove the cuff from your wrist and press the „START“ button to stop the inflation process.
- Too frequent and successive measurements may lead to disturbances of the blood circulation and cause injury.
- If you experience any discomfort during measurement, such as pain in the arm or other complaints, press the „START“ button to deflate the cuff. Release the cuff and remove it from your arm.
- Please ensure that the operation of the device does not permanently impair the patient's blood flow.
- Do not use the cuff on damaged skin.
- The material of the cuff has been tested and meet the requirements for the biological evaluation of medical devices in accordance with the standards DIN EN ISO 10993-5 and DIN EN ISO 10993-10. The constituent materials do not have the potential to produce any irritation or allergic reactions.
- If you need to test the measuring accuracy of the Geratherm® active control, please contact the manufacturer or your retailer.

## Cleaning

- Please use a soft cloth for cleaning and do not use any abrasive or volatile cleaning agents. Use a solvent-free cleaner.
- Do not wash the cuff in the washing machine or dishwasher!
- The service life of the cuff can vary according to the frequency of washing, the condition of the skin and the manner of storage. The typical service life is approximately 10.000 measurements.
- Do not clean the blood pressure monitor while it is charging. Before cleaning, always remove it from the charger first.

EN

## Maintenance and Disposal

- **Warning:** Do not carry out any maintenance or upkeep while the device is in operation.
- Dispose of the unit, accessories and components in accordance with local regulations.
- In the event of malfunction, do not attempt to repair the device yourself. Only permit repairs to be carried out by authorised technicians.
- Upon request, the manufacturer will provide circuit diagrams, parts lists, descriptions, calibration instructions, etc., to support service personnel in the repair of parts.
- The device does not need to be calibrated for two years after reliable maintenance.
- Do not dispose your blood pressure monitor in fire. The rechargeable battery can explode and lead to injury or death.
- Batteries (battery pack or inserted batteries) must not be exposed to excessive heat such as sunshine, fire or similar.

## Storage

- When it is not in use, please store the blood pressure monitor in a dry room and protect it from extreme humidity, heat, fluff, dust and direct sunlight. Do not place any heavy objects on top of the storage bag.
- This device contains sensitive parts and must be handled with care. Pay attention to the storage and operating instructions.
- It takes at least 30 minutes for the device to warm up from the minimum storage temperature to be ready for its intended use. The device must be allowed to cool down from the maximum storage temperature for at least 30 minutes before it is ready to be used again.
- Store the blood pressure monitor in a cool, dry and well-ventilated place. Keep it away from fire and heat sources, otherwise the battery could explode!

EN

## Warranty

The warranty for this blood pressure monitor is valid for any error on the part of the manufacturer under normal use for 2 years from the date of purchase.

If your unit does not function properly due to defective parts or assembly, we will repair it free of charge. With the exception of the battery and cuff, all parts of the unit are subject to this warranty. Damage caused by improper use or handling of this blood pressure monitor is not covered by the warranty.

## Information You Should Know Before Operating the Unit

### What is blood pressure?

A force is created by the heart as the ventricle forcibly ejects blood into the blood vessels and through the vascular system. Another force is created by the arteries as they resist the blood flow. Blood pressure is the result of these two forces.


### What are systolic and diastolic blood pressures?

Systolic blood pressure is the highest pressure at the heart's maximum contraction. Diastolic blood pressure is the lowest pressure when the heart is resting.

## Information You Should Know Before Operating the Unit

### Is my blood pressure normal?

See the following blood pressure classification chart released by the WHO (World Health Organization) for evaluation of your blood pressure level.

EN

| Classification of blood pressure monitor | Systolic mmHg | Diastolic mmHg | Colour indicator |
|------------------------------------------|---------------|----------------|------------------|
| Optimal                                  | < 120         | < 80           | green            |
| Normal                                   | 120 - 129     | 80 - 84        | green            |
| High - Normal                            | 130 - 139     | 85 - 89        | green            |
| Grade 1 Hypertension                     | 140 - 159     | 90 - 99        | yellow           |
| Grade 2 Hypertension                     | 160 - 179     | 100 - 109      | orange           |
| Grade 3 Hypertension                     | ≥ 180         | ≥ 110          | red              |


## Information You Should Know Before Operating the Unit

### What does low blood pressure mean?

In general, a lower blood pressure reading is better unless it causes unpleasant symptoms such as fainting and/or lightheadedness.

### Fluctuation and variation in blood pressure

The following chart shows possible blood pressure fluctuations during a 24-hour period.

The following factors will influence blood pressure measurement results and cause variations:

- Bathing
- Drinking alcohol
- Moving
- Meals
- Thoughts
- Conversation
- Exercise
- Stress
- Temperature change
- Smoking etc.

## Measurement Principle

This product uses the oscillometric measuring method to detect blood pressure. Before every measurement, the unit establishes a "zero pressure" equivalent to the air pressure.

Then it starts inflating the cuff. As it does so, the unit detects pressure oscillations generated by the pulsation of the blood flow.


The device compares the longest and the shortest time intervals of detected pulse waves with the mean time interval and then calculates the standard deviation.

Along with the blood pressure reading, the device also shows a symbol (arrhythmia icon) if the heart beat is irregular.

EN


## About the Unit

### Information about the Device

### Content:


1. Blood pressure monitor wrist type (GT-1215)
2. User manual
3. USB cable
4. Storage bag

The cuff is suitable for wrist circumferences between 13.5 and 21.5 cm.

## About the Unit

### Explanation of the Display


| Symbol       | Description              | Explanation                                                                  |
|--------------|--------------------------|------------------------------------------------------------------------------|
| <b>SYS</b>   | Systolic blood pressure  | High pressure result                                                         |
| <b>DIA</b>   | Diastolic blood pressure | Low pressure result                                                          |
| <b>PULSE</b> | Pulse                    | Pulse/Minute; Beats/ Minute                                                  |
|              | Movement detector        | Movement will lead to inaccurate readings                                    |
|              | Arrhythmia               | Irregular heartbeat Detection                                                |
|              | Position control         | Taking up the correct position is necessary for obtaining accurate readings. |
| <b>AVG</b>   | The average value        | The average value of the latest 3 blood pressure measurement results         |
|              | Low battery              | Battery is low and needs to be charged                                       |
|              | WHO Grade                | Classification of the reading according to the WHO                           |
|              | Current date and time    | Month/day, hour/minute                                                       |
|              | Memory                   | The displayed measurement values is from the memory.                         |
|              | Heartbeat                | Heartbeat detection during the measurement                                   |

EN

## Preparation for Use

### Power Supply and Charge Power

1. The battery of the Geratherm® active control is a built-in lithium-polymer battery with an electrical charge of 420 mAh.
2. Please use a power adapter with a USB connection (not included) or another source of power with a USB connection and the enclosed USB cable for charging the rechargeable battery, as shown in the following illustrations:

**⚠ Note:** **Optional Adapter**  
(Please use an authorised adapter)  
Input: AC 100 - 240 V ~ 50/60 Hz 0.2 A Max  
Output: 5 V — 1000 mA

**⚠ Power needs to be charged under the following circumstances:**

→ **Lo** on the LCD display

## Preparation for Use

If you switch on the blood pressure monitor and the LCD display flashes "  ", this means that the battery charge is low. Please recharge it in good time. You can still carry out the measurement.

If the LCD displays "  L  " means the battery is too low, the blood pressure monitor will turn off automatically, you must charge the power at once.

During the process of charging, the LCD display will display the blinking power level, just like the following pictures:

When there is no power, the LCD will displays as the following when it is charging:

When there is one level of power, the LCD will displays as the following when it is charging:

When there are two levels of power, the LCD will displays as the following when it is charging:

When there are three levels of power, it means the power is full, the LCD will displays as the following when it is charging:

EN

## Preparation for Use

### CAUTION:

- The battery of Geratherm® active control is a built-in rechargeable lithium-polymer battery, please do not attempt to disassemble the blood pressure monitor or force open the built-in battery by the unauthorized maintenance personnel.
- In normal use, the battery can be recharged about 300 times. If the battery fails to recharge or the unit cannot be used normally, please contact your dealer. If you measure your blood pressure three times a day, the unit can be used for up to 15 days without recharging.
- Do not attempt to replace the battery of your blood pressure monitor. The battery is built-in and is not replaceable
- Avoid recharging your blood pressure monitor in extremely high or low ambient temperatures.
- Do not clean the blood pressure monitor during recharging. Always remove the charging unit before cleaning.
- While it is charging, do not touch the charger connection and the patient simultaneously.

EN

## Preparation for Use

### Activate your blood pressure monitor

Your blood pressure monitor is activated when the date and time are set.

### Setting Date and Time

It is important to set the clock before using your blood pressure monitor, so that a time stamp can be assigned to each record that is stored in the memory. (time format: 24 h)

1. When the monitor is off, hold pressing "MEM" button about 3 seconds to enter the mode for year setting.

2. Press "START" button to change the year. Each press will increase the numeral by one in a cycling manner.


3. When you get the right year, press "MEM" button to set down and turn to next step.


EN

## Preparation for Use


4. Repeat steps 2 and 3 to set the [MONTH] and [DAY].

5. Repeat steps 2 and 3 to set the [HOUR] and [MINUTE].

6. Repeat steps 2 and 3 to switch position control on or off. After confirming this as described in step 3, all the settings will be successively shown on the display. Then the unit will turn off.

## Preparation for Use

### Attaching the cuff

1. Wrap the cuff around the bare wrist. The display should face the side of the palm of the hand.
2. Fasten the cuff. Make sure it is not too tight. The cuff's upper edge should be approximately 1 cm (0.39 inch) from the palm line.

EN

### Measuring posture

1. Rest for five minutes before measuring. Wait at least three minutes between measurements. This will enable your circulation to recover.
2. Sit upright and position your lower arm in such a way that the unit is at heart level. Relax and retain a natural posture during measurement.
3. Measure and record blood pressure at the same time every day to establish your blood pressure pattern.


## How to Operate the Unit


### Measuring Blood Pressure

1. When the monitor is off, press the "START" button to turn on the monitor, and it will finish the whole measurement.

LCD display

0 – value is shown

Inflating and  
measuring

Results are displayed  
and stored  
in the memory


2. As soon as the measurement is completed, the blood pressure and pulse readings appear on the display. Press the "START" button to switch off the unit. Otherwise, the unit will switch off automatically after one minute.

## How to Operate the Unit

### **Note:**


If you have turned on the Position Control function, once the measurement starts, the blood pressure monitor will detect the gesture first. The wrist must be between the angle of 30° and 45°. If the position is out of this angle, the measurement will not start and the display will show Symbol  + ERR until you have the correct position.

EN

### **Note:**

If you have the correct position, the display will show  for about 3 seconds, and then continue with the measurement.

## How to Operate the Unit

### Retrieving Readings

1. When the monitor is off, please press "MEM" button to show the average value of the latest three records.

(Note: If the records are less than 3 groups, the LCD will display the recent record instead.)

EN


**Note:**

The device can record 60 measurements.

2. By pressing again on the "MEM" button, all stored readings can be retrieved. The number of the stored reading is shown in the bottom left-hand corner. This appears alternately with the time when the reading was taken. The date is shown underneath.

The most recent record (1) is shown first. Each new measurement is assigned to the first (1) record. All other records are pushed back one digit (e.g., 2 becomes 3, and so on), and the last record is dropped from the list.

## How to Operate the Unit

### Deleting Readings

1. Hold pressing "MEM" button for 3 seconds when the monitor is in memory mode, the display will show "dEL ALL" ("delete all") .

2. Press "MEM" button to confirm deleting and the monitor will display "dEL dOnE" ("deletion done") and then turn off.


3. If there is no record press "MEM" button, the LCD will display "0".

#### **Note:**

To exit out of delete mode without deleting any records, press "START" button before pressing "MEM" button to confirm any delete commands.

EN

|                                                                                                                  |                                                                                   |
|------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|
| <p>Do not drop the unit.</p>                                                                                     |  |
| <p>Do not modify or disassemble the unit or the cuff.</p>                                                        |  |
| <p>Do not twist the cuff.</p>                                                                                    |  |
| <p>Use a cloth moistened with water or neutral detergent to clean the body of the unit and then wipe it dry.</p> |  |

## Care and Maintenance

Avoid thinner, benzine, and other harsh cleaners.

Keep the unit in a suitable place.  
Avoid high temperature, direct sunlight, high moisture, and dust.

Do not press the "START" button until the cuff has been properly fastened on the wrist.

EN

## Error messages

EN

| Problem           | Display symbol                                                                              | Cause                                                                 | Solution                                                                                                                    |
|-------------------|---------------------------------------------------------------------------------------------|-----------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| No power          | Nothing appears on the display                                                              | Battery charge too low                                                | Charge battery                                                                                                              |
| Low battery power |  Indicator | Battery charge low                                                    | Charge battery                                                                                                              |
| No power          | Nothing appears on the display                                                              | Adapter / charging cable is not connected correctly                   | Check adapter and charging cable and plug in correctly                                                                      |
| No display        | No display when buttons pressed                                                             | Unit is not activated                                                 | Press and hold down „MEM“ button to activate the unit                                                                       |
| Error messages    | E 1                                                                                         | The cuff is not properly fastened or pressure in the cuff is too high | Rest for a short period of time, re-fasten the cuff and measure the blood pressure again                                    |
|                   | E 2                                                                                         | The unit has detected movement or pulse rate is too low               | Movement can influence the measurement. Relax for a moment and then measure again                                           |
|                   | E 3                                                                                         | No pulse signal is detected                                           | Loosen clothing around the arm and measure again                                                                            |
|                   | E 4                                                                                         | The measuring process was not successful                              | Rest for a short time and retake the measurement                                                                            |
|                   | EExx                                                                                        | A calibration error has occurred                                      | Repeat the measurement. If the problem persists, contact your dealer or our customer service department for further support |
|                   |            | Measured values are outside the measurement range                     | Relax for a moment. Re-fasten the cuff and measure again. If the problem re-occurs, consult your physician                  |

## Specifications

|                                    |                                                                                                                                                                                                                                            |
|------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Model no.                          | GT-1215                                                                                                                                                                                                                                    |
| Power supply                       | 3.7 V 420 mAh Built-in rechargeable lithium-polymer battery                                                                                                                                                                                |
| Display mode                       | Digital LCD Display V.A. 44.8 mm x 35.6 mm                                                                                                                                                                                                 |
| Measurement mode                   | Oscillographic testing mode                                                                                                                                                                                                                |
| Measurement range                  | Rated cuff pressure: 0mm Hg ~ 299 mmHg (0 kPa ~ 40 kPa)<br>Measurement pressure:<br>SYS: 60 mmHg ~ 230 mmHg (8.0 kPa ~ 30.7 kPa)<br>DIA: 40 mmHg ~ 130 mmHg (5.3 kPa ~ 17.3 kPa)<br>Pulse value: (40 - 199) beat/minute                    |
| Accuracy                           | Pressure: $\pm 3$ mmHg (0.4 kPa)<br>Pulse value: $\pm 5$ %                                                                                                                                                                                 |
| Operating condition                | A temperature range between: +5 °C and +40 °C<br>A relative air humidity range of: $\leq 85$ % non-condensing but not requiring water vapour partial pressure of more than 50 hPa<br>An atmospheric pressure range of: 700 hPa to 1060 hPa |
| Storage & transportation condition | Temperature: -20 °C to 60 °C<br>A relative air humidity range of $\leq 93$ %, non-condensing, with water vapour pressure up to 50 hPa                                                                                                      |
| Cuff size                          | 13.5 cm ~ 21.5 cm                                                                                                                                                                                                                          |
| Weight                             | 104 g (including the cuff)                                                                                                                                                                                                                 |
| External dimensions                | approx. 84 mm x 70 mm x 40 mm                                                                                                                                                                                                              |
| Attachment                         | USB cable, storage bag, user manual                                                                                                                                                                                                        |
| Degree of protection               | IP22, The device is protected against penetration by solid objects $\geq 12.5$ mm in diameter and against the ingress of water in the form of vertically falling drops.                                                                    |
| Device classification              | Battery mode: ME device with built-in lithium-polymer rechargeable battery                                                                                                                                                                 |
| Software Version                   | A01                                                                                                                                                                                                                                        |

EN

Subject to changes in the interests of technical progress.

## Quality Standard

Geratherm® is certified in accordance with Council Directive 93/42/EEC and EN ISO 13485 and is entitled to affix the CE-mark C€0197 (Notified Body: TÜV Rheinland LGA Products GmbH).

## Symbol Index

|    |                                                                                   |                                                         |                                                                                   |                                                             |
|----|-----------------------------------------------------------------------------------|---------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------------------------|
| EN |  | Follow the instructions for use                         |  | Type BF Equipment                                           |
|    |  | Keep dry                                                |  | Manufacturer                                                |
|    |  | Store at relative humidity levels ≤ 93 %                |  | Batch code<br>(YYMMXXX, Year/Month / Identification number) |
|    |  | Store between -20 °C and +60 °C                         |  | Caution! Read the instruction manual.                       |
|    |  | The device must not be disposed of with household waste |  | Manufacturing date                                          |
|    |  | Direct current                                          |  | Reference number                                            |

## Information on electromagnetic compatibility (EMC)

Electronic devices such as PCs and mobile phones can lead to medical devices in operation being exposed to electromagnetic interference from other devices. This can lead to malfunction of the medical device and create a potentially unsafe situation.

Medical devices should also not interfere with any other equipment.

The EN 60601-1-2 standard regulates the requirements for EMC (electromagnetic compatibility) and defines the levels of immunity to electromagnetic interference and the maximum electromagnetic emission levels for medical devices.

This blood pressure monitor, which is manufactured by Geratherm Medical AG, complies with the EN 60601-1-2 standard in relation to both immunity and emissions.

Nevertheless, special precautions should be observed:

Please only use the device indoors and not in the vicinity of mobile phones or microwave ovens. This device should not be used adjacent to or stacked with other devices; and if it has to be used adjacent to or stacked with other devices, it should be observed to verify normal operation in the configuration in which it is being used.

In the case of devices exceeding the maximum power of 2 W, the minimum distance from your blood pressure monitor should be 3.3 metres.

### Warning:

The use of accessories, transducers and cables other than those specified or provided by the manufacturer of this device may result in increased electromagnetic emissions or reduced electromagnetic immunity of the device and may cause it to malfunction.

EN

| Guidance and manufacturer's declaration – electromagnetic emissions                                                            |  |            |                                                                                                                                                                                                                          |
|--------------------------------------------------------------------------------------------------------------------------------|--|------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments: |  |            |                                                                                                                                                                                                                          |
| Emissions test                                                                                                                 |  | Compliance | Electromagnetic environment - guidance                                                                                                                                                                                   |
| RF emissions CISPR 11                                                                                                          |  | Group 1    | RF energy is used only to maintain device's operation. Therefore, its RF emissions are so low that it's not likely to cause any interference in nearby electronic equipment.                                             |
| RF emissions CISPR 11                                                                                                          |  | Class B    | The device is suitable for use in all establishments, including domestic establishments, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2                                                                                               |  | Class A    |                                                                                                                                                                                                                          |
| Voltage fluctuations / flicker emissions IEC 61000-3-3                                                                         |  | complies   |                                                                                                                                                                                                                          |

| Guidance and manufacturer's declaration – electromagnetic emissions                                                                                                                                                           |                              |                              |                                                                                                                                                |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------|------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| The device or system is intended for use in the electromagnetic environments listed below. The customer and/or user of the device or system must ensure that it is used in an electromagnetic environment as described below: |                              |                              |                                                                                                                                                |
| Immunity test                                                                                                                                                                                                                 | IEC 60601 – test level       | Compliance level             | Electromagnetic environment - guidance                                                                                                         |
| Electrostatic discharge (ESD) IEC 61000-4-2                                                                                                                                                                                   | ± 6 kV contact<br>± 8 kV air | ± 6 kV contact<br>± 8 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8                                                                                                                                                                       | 3 A/m                        | 3 A/m                        | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.      |

| Recommended separation distances between portable and mobile RF communication                                                                                                                                                                                                                                                                       |                                         |                                         |                                          |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------|-----------------------------------------|------------------------------------------|
| The device is intended for use in an electromagnetic environment where radiated RF disturbances are under control. User can help prevent electromagnetic interference by keeping the device at a minimum distance from portable and mobile RF communications equipment (transmitters). Below table details the maximum output power of transmitter: |                                         |                                         |                                          |
| Rated maximum output power of transmitters in Watt                                                                                                                                                                                                                                                                                                  | Separation distance / m                 |                                         |                                          |
|                                                                                                                                                                                                                                                                                                                                                     | 150 kHz to 80 MHz<br>$d = 1.2 \sqrt{P}$ | 80 MHz to 800 MHz<br>$d = 1.2 \sqrt{P}$ | 800 MHz to 2.5 GHz<br>$d = 2.3 \sqrt{P}$ |
| 0.01                                                                                                                                                                                                                                                                                                                                                | 0.12                                    | 0.12                                    | 0.23                                     |
| 0.1                                                                                                                                                                                                                                                                                                                                                 | 0.38                                    | 0.38                                    | 0.73                                     |
| 1                                                                                                                                                                                                                                                                                                                                                   | 1.2                                     | 1.2                                     | 2.3                                      |
| 10                                                                                                                                                                                                                                                                                                                                                  | 3.8                                     | 3.8                                     | 7.3                                      |
| 100                                                                                                                                                                                                                                                                                                                                                 | 12                                      | 12                                      | 23                                       |
| For transmitters rated at a maximum output power not listed above, the recommended separation distance $d$ in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.                |                                         |                                         |                                          |
| NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.                                                                                                                                                                                                                                                      |                                         |                                         |                                          |
| NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.                                                                                                                                                                                 |                                         |                                         |                                          |

| Guidance and manufacturer's declaration – electromagnetic immunity                                                                                                                                                            |                             |                  |                                                                                                                                                                                                                                                                                                                                                                              |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------|------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The device or system is intended for use in the electromagnetic environments listed below. The customer and/or user of the device or system must ensure that it is used in an electromagnetic environment as described below: |                             |                  |                                                                                                                                                                                                                                                                                                                                                                              |
| Immunity test                                                                                                                                                                                                                 | IEC 60601 – test level      | Compliance level | Electromagnetic environment - guidance                                                                                                                                                                                                                                                                                                                                       |
|                                                                                                                                                                                                                               |                             |                  | Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.                                                                                                                                  |
|                                                                                                                                                                                                                               |                             |                  | Recommended separation distance:                                                                                                                                                                                                                                                                                                                                             |
| Conducted RF IEC 61000-4-6                                                                                                                                                                                                    | 3 Vrms<br>150 kHz to 80 MHz | 3 Vrms           | $d = 1.2 \sqrt{P}$                                                                                                                                                                                                                                                                                                                                                           |
| Radiated RF IEC 61000-4-3                                                                                                                                                                                                     | 3 V/m<br>80 MHz to 2.5 GHz  | 3 V/m            | $d = 1.2 \sqrt{P}$<br>80 MHz to 800 MHz                                                                                                                                                                                                                                                                                                                                      |
|                                                                                                                                                                                                                               |                             |                  | $d = 2.3 \sqrt{P}$<br>800 MHz to 2.5 GHz                                                                                                                                                                                                                                                                                                                                     |
|                                                                                                                                                                                                                               |                             |                  | where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).<br>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup> |
|                                                                                                                                                                                                                               |                             |                  | Interference may occur in the vicinity of equipment marked with the following symbol:                                                                                                                                                                                                                                                                                        |
|                                                                                                                                                                                                                               |                             |                  |                                                                                                                                                                                                                                                                                             |
| NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.                                                                                                                                                            |                             |                  |                                                                                                                                                                                                                                                                                                                                                                              |
| NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.                                                           |                             |                  |                                                                                                                                                                                                                                                                                                                                                                              |

| Guidance and manufacturer's declaration – electromagnetic immunity |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|--------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| a)                                                                 | Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device. |
| b)                                                                 | Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| 1)                                                                 | This equipment needs to be installed and put into service in accordance with the information provided in the user manual;                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| 2)                                                                 | Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance $d = 3,3$ m away from the equipment. (Note: As indicated in Table 6 of IEC 60601-1-2 for ME EQUIPMENT, a typical cell phone with a maximum output power of 2 W yields $d = 3,3$ m at an IMMUNITY LEVEL of 3 V/m)                                                                                                                                                                                                                                                   |

The current version of the standards is valid.

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|                                                                                                                                   |                                                                                                                                                                                                                                                                                                                  |
|-----------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>Non far cadere il dispositivo.</p>                                                                                             |  An illustration of the ConTherm device falling from a height, indicated by motion lines. A large prohibition sign (a circle with a diagonal line) is overlaid on the image to indicate that dropping the device is prohibited. |
| <p>Non modificare o smontare il dispositivo o il polsino.</p>                                                                     |  An illustration showing a screwdriver being used to attempt to open or modify the device. A large prohibition sign is overlaid on the image to indicate that disassembly or modification is prohibited.                        |
| <p>Non torcere il polsino.</p>                                                                                                    |  An illustration of the device being twisted or bent. A large prohibition sign is overlaid on the image to indicate that twisting the device is prohibited.                                                                     |
| <p>Per pulire la custodia utilizzare un panno umidificato con acqua o con detergente neutro e successivamente asciugare bene.</p> |  An illustration of a hand using a cloth to clean the device. This illustrates the correct cleaning procedure.                                                                                                                  |

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