

ZEUS-200S/150



Output Specification

Mode	ZEUS-200S	ZEUS-150
Pure Cut	200W / 200Ω	150W / 200Ω
Blend Cut 1	150W / 200Ω	120W / 200Ω
Blend Cut 2	120W / 200Ω	100W / 200Ω
Contact Coagulation	80W / 200Ω	80W / 200Ω
Spray Coagulation	80W / 500Ω	60W / 500Ω
Bipolar Coagulation	80W / 100Ω	80W / 100Ω

- Main Frequency : 400, 500KHz
- Power Source : AC120V or AC230V, 50Hz / 60Hz
- Protection Class : Class 1, Type CF
- Dimension(WxDxHmm) : 330×320×150
- Weight : 8.0Kg

Feature

• Various Function with Compact Design

Cut(Pure Cut, Blend 1, Blend 2), Coag.(Contact, Spray), Bipolar Coag.; Fully enough for most of the operations.

• Best Safety with Feedback System

REM(Return Electrode Monitoring) function monitors size of contact area between a patient and the pad. If contact is inappropriate, it automatically blocks the high-frequency current to minimize the danger of burning incident.

• Different Tone Audible Sound

Each working condition can be checked by sound and indication lamp
The level if alarm sound can be adjusted during an operation of application mode.

• Auto Select Function

Without pressing cut or coag. button, use twin button handle or double foot switch for working cut or coag. mode.
Press each mode button for 3 sec., output power will be set 1W, it enables to adjust from 1W to maximum power quickly.

Standard Accessories

Code No.	Name of Item				
FS02-01R	Double Foot Switch				
HD02-01D	Disposable Twin Button Handle				
HD01-01R	Monopolar Handle & Cable				
PL02-04R	Silicon Patient Plate(240X150mm)				
PL03-02D	Disposable Patient Plate (2ea)				
CA02-03R	Patient Return Plate Cable				
CA01-01R	Bipolar Forceps Cable				
BF01-01R	Bipolar Forceps(Bayonet 17.8cm 1.0 Tip)				
EL01-02D	Electrode(Knife 2.4*70mm)				
EL02-02D	Electrode(Ball 5mm)				
EL03-02D	Electrode(Loop 6mm)				
EL04-02D	Electrode(Needle 2.4*70mm)				
EL04-04D	Electrode(Needle Angled 2.4*70mm)				
		EL01-02D / EL02-02D EL03-02D / EL04-02D / EL04-04D			