

70L CAPACITY

YELLOW PLASTIC BIO MEDICAL WASTE BIN

Safe. Hygienic. Responsible.

Designed for safe collection and disposal of biomedical waste in hospitals, clinics, laboratories, and healthcare facilities.



FOOT PEDAL OPERATION

Hands-free operation ensures maximum hygiene and safety.



DURABLE CONSTRUCTION

Made from high-quality HDPE plastic, impact resistant and long-lasting.



LARGE CAPACITY

70L capacity suitable for high volume waste disposal.



EASY TO CLEAN

Smooth surface for easy cleaning and disinfection.



LEAK PROOF

Sealed-design prevents leakage



CORROSION RESISTANT

Resistant to chemicals and moisture



ENVIRONMENT FRIENDLY

Made from recyclable materials

PRODUCT VIEWS

FRONT VIEW



BACK VIEW



OPEN VIEW



PEDAL DETAIL



SPECIFICATIONS

Product Name	Yellow Plastic Bio Medical Waste Bin
Capacity	70 Liters
Dimensions (L x W x H)	65 x 50 x 42 cm
Material	High Quality HDPE
Weight	4 kg
Operation	Foot Pedal
Color	Yellow
Usage	Indoor / Semi-Outdoor
Ideal For	Hospitals, Clinics, Labs, Nursing Homes
Features	Durable, Hygienic, Easy to Clean

IDEAL FOR



Hospitals



Clinics



Laboratories



Nursing Homes



Healthcare Centers



SAFE DISPOSAL



HYGIENIC DESIGN



EASY MAINTENANCE



SUSTAINABLE CHOICE



MEDISPOT
MEDICAL SUPPLIES

Trusted Partner in Healthcare Solutions

Providing high-quality medical products for a safer and healthier environment.

✉ info@medispotmedical.com
🌐 www.medispotmedical.com